

SOCOREX® *ultra* 1810/2810

ANIMAL HEALTH SYRINGE VETERINÄR SPRITZEN SERINGUES VÉTÉRINAIRES

OPERATING INSTRUCTIONS BEDIENUNGSANLEITUNG MODE D'EMPLOI

! SAFETY PRECAUTIONS

- Read the operating instructions carefully before use. Keep available for future reference
- Observe the manufacturer's recommendations
- Follow instructions for proper use of injectables
- Before each use, check instrument for perfect working conditions
- Purge, rinse, dry syringe and lubricate the tightness parts carefully after each use
- Never try to disassemble blocked parts and do not use force to separate plunger from cylinder or valves
- Repeated use of the same needle may carry over contaminations to other animals
- Instrument manufacturer cannot be made responsible for needle breakage
- Syringe should be used by professionals only
- Mind possible hand fatigue and its medical consequences during prolonged vaccine campaign
- Do not use syringe for purposes other than intended for
- A saline environment may hasten replacement of instruments parts
- Do not force on handle if needle is blocked. High pressure inside syringe, may lead to glass cylinder breakage and valve damages
- Only use original spare parts

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DESCRIPTION

Congratulations on buying a high quality Swiss product. Each syringe from the *ultra* family allows trouble free injection of various water or oil based vaccines with excellent reproducibility. The syringe with spring loaded plunger and three-way valve system is intended for micro-range injections (up to 0.1, 0.2, 0.3 and 0.5 mL) such as for fishes, chicks and small laboratory animals. Larger volumes (1 and 2 mL) are ideal for poultry or pigs.

It can be autoclaved fully assembled at 121°C (20 minutes). Untighten valves and casing nut prior autoclaving.

Socorex® *ultra* syringe is available in 5 models:

- Model 1810.2. with feed tubing,
- Model 1810.1. with vial holder,
- Model 1810.6. with feed tubing and Abacus™ shot counter,
- Twin model 2810.2 with feed tubing,
- Twin model 2810.1 with vial holder.

Each feed tubing syringe can be transformed into a vial holder version, and vice versa, with a respective conversion set (available separately as accessory).

Valves are self-contained, avoiding loss of small parts. Volume is easily set with a micrometer screw and locked by nut. The graduated glass barrel is interchangeable.

Veterinary needles and drenchers to be ordered separately. Details on Luer Lock stainless steel needles available upon request.

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HOLDING THE SYRINGE

Ergonomic, user friendly handle, perfect balance and smooth plunger movement, ensure fatigue-free operation even during mass vaccination at rates of up to 1000 injections per hour.

Place hand on handle and move plunger with hand palm.

To change position of syringe casing (12), untighten screw (7), pull toward the body, rotate the syringe body into one of the eight possible slots. Push back body and retighten screw (7) keeping plunger in forward position.

OPERATION

The instrument is supplied ready to use, no assembling required.

Self-refilling syringe with feed tubing

- Aspiration from bottle at belt:
Puncture bottle cap with feed needle and vent cannula and connect feed tubing to the feed needle. Filling by gravity the entire silicone tubing, facilitate priming. Verify absence of rubber residue in tubing before connection. Repeated puncture of the rubber seal is not recommended.
- Aspiration from bottle hanging from holder (neck toward the bottom):
Connect feed tubing to the optional draw-off needle with air inlet and filter (ref. 1.187.590).
- Aspiration from open bottle:
Connect sinker to the end of the feed tubing and immerse into liquid.

Self-refilling syringe with vial feeding

To draw liquid out of an original vaccine bottle with septum closure, follow below steps :

1. Remove bottle lock (Fig. 1).
2. Insert bottle into lock (Fig. 2, max $\varnothing$ 20 mm).
3. Holding syringe upside down, pierce septum with needle. Screw the lock back on (Fig. 3). Do not over-tight- risk of valve cartridge damage.
4. Turn syringe and activate it.

Maximum bottle capacity on the Twin syringe is either 2 x 50 mL bottles or 1 x 100 mL and 50 mL bottle.


Fig. 1


Fig. 2


Fig. 3


Fig. 4

Priming

Set instrument on maximum volume. Prime feed tubing by activating plunger, until liquid completely fills the syringe. Hold syringe upright to expel all air bubbles and/or residues by activating plunger to minimum of five times. The lost liquid during priming can be collected in a clean recipient and returned to container. After changing bottle, verify proper priming (no air bubble) to guarantee optimal performances.

Volume setting (Fig 5 and 6)

Activate plunger a few times before setting volume to enable free plunger movement (stick slip).

Setting volume on models with capacity of up to 0.3 mL

Loosen locking nut (5) and turn socket (4) to center O-ring on desired graduation on models 0.1, 0.2 and 0.3 mL (Fig 5).

Setting volume on models with capacity of 0.5 - 2 mL

For 0.5, 1 and 2 mL, center plunger edge on desired graduation (Fig 6).

Tighten locking nut (5) in order to prevent any volume alteration. Activate a few times the plunger to verify good operation.

On 2810 Twin syringe, the volume is independently adjustable on each barrel.

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
Injection

Attach injection needle (to be ordered separately) to nozzle (mind adequate needle size). During injection, do not apply force on plunger rod and check output liquid from needle.

Too much pressure inside barrel can cause barrel breakage or leakage (risk of injury).

Optional manifold (Cat. No. 1.287.685) for Twin syringes allows injection of two different liquids through the same needle.

Manifold mounting (Twin model 2810)

Loosen both casing nuts (7). Position the manifold onto both Luer lock nozzles (16). Tighten both manifold security rings. Tighten both casing nuts.

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When dosing is completed, remove feed tubing and activate plunger several times to eliminate any residual liquid. Rinse, disassemble, clean and dry the instrument parts according to the instructions in the "Maintenance" chapter.

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To remove the bottle from vial feeding model, turn syringe upside down and unscrew the lock (Fig. 4). Activate plunger to remove residual liquid.

Never store a syringe without prior purging, cleaning and drying parts. Lubrication of tightness parts (plunger O-ring, barrel and valves) avoids blockage over time. Medical grade white and silicone oils recommended. Other stable food oils can also be used.

STERILIZATION

The *ultra* syringe can be autoclaved fully assembled at 121°C, for 20 minutes in a sterilization pouch (Qualitix cat. No. 3210.1326) laying directly on a metal grid.

Do not sterilize an instrument that has not been previously purged and cleaned. See chapter "Maintenance".

Lubricate plunger O-ring, cylinder and valves before autoclaving the syringe to avoid blockage. See chapter "End of dosing".

Slightly untighten both valve connectors (15 and 16) and casing nut (7).

After sterilization, the instrument must be cool and completely dry prior to retightening all connectors and nuts. Correct autoclaving procedure is the sole user responsibility.

Abacus™ shot counter is not autoclavable and must be removed from the syringe prior to sterilisation.

Alternatively, syringe can be sterilised in boiling water at 100 °C for 30 minutes.

MATERIALS

All parts in contact with liquid are made of carefully selected and chemically inert materials:

- Feed tubing: silicone
- Barrel: sodocalcic AR glass (0.5, 1, 2 mL) and borosilicate glass 3.3 (0.1 to 0.3 mL)
- Plunger: stainless steel
- Barrel washer, plunger O-ring: FPM (FDA compliant)
- Valve: PP, stainless steel and FPM (FDA compliant)
- Syringe casing: PA
- Sink and feed cannula: nickel/chrom plated brass
- Vial holder body: PA
- Needle assembly: PA, stainless steel, FPM, silicone

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MAINTENANCE

Disassembling

Remove needle and detach feed tubing or vial holder. Hold handle (6) and unscrew casing nut (7) while holding syringe body (fingers on barrel).

Slowly pull forward casing to disengage plunger (8). Never apply excessive force if the plunger is blocked. Remove barrel, watch out not loosing barrel washer (11) Untighten both connectors (15 and 16), remove valves (13 and 14)

Cleaning

Syringe and accessories must be carefully cleaned after each use.

Aspirate water or an appropriate cleaning solution (detergent) and rinse thoroughly. Solution should not contain any residue.

Purge cleaning solution by activating plunger several times. Use isopropanol only, avoid using ethanol or methanol. It is mandatory to let dry all parts (valves, connectors, syringe casing, plunger and barrel) before reassembling.

Valves reassembling (Fig 7)

Introduce each valve in its specific housing, respecting mechanical coding. Valves for syringes with volumes up to 0.5 mL are red; for volumes above 0.5 mL, they are yellow. Lubricate valve if syringe is stored.

Screw inlet and outlet valves. Do not use tool, neither apply excessive force while manually tightening.

Plunger O-ring replacement (Fig. 8 and 9)

Applies only to 0.1, 0.2 and 0.3 mL models. Use thin needle to remove O-ring. Avoid damaging the O-ring groove.

Mount new lubricated O-ring on guide (1). Push O-ring into the plunger groove (2).

Plunger reassembly (Fig. 10 and 11)

Slightly lubricate barrel entrance (1), insert half of plunger into barrel without force (2). Place barrel into the syringe casing. Couple plunger onto plunger rod (3) and push the plunger fully inside barrel (4).

Introduce the syringe casing into the desired handle slot and partially screw casing nut (7). Position the barrel as to see the graduation. Push the plunger fully into the barrel by pressing on the piston rod (1) and tighten the nut (7) without excessive force. Assure correct syringe operation with several plunger movements. Plunger should move freely and clearly stop at end of travel.


Fig. 7

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
Fig. 11


Fig. 12

TROUBLESHOOTING	POSSIBLE CAUSE	REMEDY
Erratic operation	Syringue used without liquid inside Plunger, O-ring and barrel not properly lubricated	Use with liquid Use appropriate lubricant
Blocked valve	Syringe not cleaned Valve parts blocked by residue	Disassemble and clean Soak valve and clean (remove residue)
Leaking valve	Residue in valve (i.e. rubber residue from septum, cristalization) prevent tightness Damaged valve	Clean valve Replace part
Wrong volume Inadequate performance	See leaking valve Priming not completed, air bubble in system Barrel washer not positioned correctly or missing	See leaking valve Prime properly Change barrel washer
Liquid accumulation behind plunger	0.1 mL – 0.3 mL models, damaged O-ring or missing 0.5, 1 and 2 mL models, this may occur after long injection series	Check plunger, change O-ring Remove liquid when doing maintenance
Plunger blocked	Syringe stored without cleaning	Do not apply force to unblock. Disassemble and soak syringe in warm water and separate plunger parts Change needle
No priming possible	Needle clogged Check for blocked valve Check for leaking valve Check for clogged needle	Clean and reassemble Change valve if required Change needle if required
Difficulty in attaching the bottle onto vial holder.	Bent needle.	Push clip and pull out needle assembly. Insert new needle assembly (Fig. 12, Cat. No 1.1810.386). ⚠ Sharp needle - risk of injuries! Do not pierce the same bottle several times.

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PERFORMANCE AND ORDERING INFORMATION

Package contents

Each *ultra* model is supplied with:

- Quick Guide
- Luer Lock nozzle
- Set of spare parts

Spare part set for syringes up to 0.3 mL includes:
2 compact cartridges (inlet valve, outlet valve), 2 barrel washers, 2 plunger O-rings, 1 guide for O-ring.

Spare part set for syringes ≥0.5 mL includes:
2 compact cartridges (inlet valve, outlet valve), 2 barrel washers.

For Twin models, the spare part set is doubled.

Additionally, vial feeding model 1810.1. and 2810.1.XXXX is supplied with:

- Vial holder
- Bottle lock
- Connection tubing

Tube feeding models, 1810.2.XXXX, 1810.6.XXXX and 2810.2.XXXX are supplied with:

- 1 m silicone tube (17)
- Sinkers (18)
- Feed needle (19) and vent cannula (20)
- Abacus™ shot counter (only for 1810.6.XXXX)

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ZUBEHÖR

Cat. No.	Description	Packaging
Valves		
1810.915	Intake valve (0.1 - 0.5 mL)	1 pc
1810.945	Outlet valve (0.1 - 0.5 mL)	1 pc
1810.916	Intake valve (1 - 2 mL)	1 pc
1810.946	Outlet valve (1 - 2 mL)	1 pc
Spare part kit		
1.1810.72001	0.1 mL	1 /pk
1.1810.72002	0.2 mL	1 /pk
1.1810.72003	0.3 mL	1 /pk
1.1810.72005	0.5 mL	1 /pk
1.1810.7201	1 mL	1 /pk
1.1810.7202	2 mL	1 /pk
Various syringe accessories		
1.187.590	Draw off needle	6 /pk
1.287.685	Single needle manifold	1 /pk
187.592	Feed needle	1 /pk
187.933	Vent cannula	10 /pk
Drenching cannulas		
376.0925S	Straight, 25 mm	3 /pk
376.1251S	Straight, 51 mm	3 /pk
376.16102S	Straight, 102 mm	2 /pk
376.0925B	Bent, 25 mm	3 /pk
376.1251B	Bent, 51 mm	3 /pk
376.16102B	Bent, 102 mm	2 /pk
Abacus™ shot counter (only for 1810 syringes)		
1.1810.962	For 0.1 / 0.2 mL syringes	1 /pk
1.1810.966	For 0.3 - 2 mL syringes	1 /pk
1.1810.160	Magnet for reset	1 /pk

ACCESSOIRES

Cat. No.	Description	Packaging
Premium lubricant ABF		
1.170.705	70 mL	1 pc
1.170.906	70 mL (pack of 6)	6 /pk
1.170.755	500 mL	1 pc
1.170.956	500 mL (pack of 6)	6 /pk
1.1810.700	White oil	1 /pc
Superior quality needles, Luer Lock *		
371.0503	0.5 x 3 mm / 25G x 1/8"	12/pk
371.0504	0.5 x 4 mm / 25G x 1/6"	12/pk
371.0603	0.6 x 3 mm / 23G x 1/8"	12/pk
371.0605	0.6 x 5 mm / 23G x 3/16"	12/pk
371.0608	0.6 x 8 mm / 23G x 1/3"	12/pk
371.0704	0.7 x 4 mm / 22G x 1/6"	12/pk
371.0706	0.7 x 6 mm / 22G x 1/4"	12/pk
371.0710	0.7 x 10 mm / 22G x 3/8"	12/pk
371.0806	0.8 x 6 mm / 21G x 1/4"	12/pk
371.0810	0.8 x 10 mm / 21G x 3/8"	12/pk
371.1010	1.0 x 10 mm / 19G x 3/8"	12/pk
371.1213	1.2 x 13 mm / 18G x 1/2"	12/pk
371.1420	1.4 x 20 mm / 17G x 3/4"	12/pk
371.1625	1.6 x 25 mm / 16G x 1"	12/pk
FishGuide™ (only for 1810 syringes)		
1.1810.015	15 mm long	1 /pk
1.1810.025	25 mm long	1 /pk
1.1810.045	45 mm long	1 /pk

* Additional sizes, see www.socorex.com

Viscosity data

Volume mL	Maximal cinematic viscosity * (cst)	Example of liquids
0.1 to 0.5	200	Glycerol 90%
1	1200	Silicon oil
2	400	Glycerol 94%
Test conditions to determine maximum viscosity		
Syringe with horizontal 1 meter feed tubing. Syringe set at maximal volume. Aspiration time < 3 sec. Needle diameter 0.8mm (21G).		

* Cinematic viscosity in centistokes (cst). Dynamic viscosity in centipoises (cps) obtained by multiplying cinematic viscosity by liquid density (kg/dm3). Performance data and viscosity limit are related to selected needle size.

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FISH INJECTIONS (ONLY FOR 1810 SYRINGES)

FishGuide™ is a convenient accessory providing for correct needle positioning while injecting fishes. It is not compatible with Twin model. The accessory is supplied in three sizes: 15 mm, 25 mm and 45 mm and does not require any tool to be installed and removed. It can turn over 360° on the syringe nozzle. The two large copper grips facilitate fish positioning. Their angle can be enlarged according to fish size.

Notes:

Improper FishGuide™ positioning or overtightening screw may damage nozzle and impair needle fitting.

Reduce risk of injuries by removing FishGuide™ when changing needle.

Assembling (Fig. 13 - 17)

1. Unscrew Luer Lock nozzle (Fig. 13).
2. Insert black jacket at a small angle ① and push until abutment ② (Fig. 14)
3. Screw Luer Lock nozzle, attach needle (Fig. 15)
4. Insert FishGuide™ ①, adjust position ②, softly tighten thumb screw ③ (Fig. 16)
5. Adjust angle as necessary (Fig. 17).


Fig. 13


Fig. 14


Fig. 15


Fig. 16


Fig. 17

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ABACUS™ SHOT COUNTER (ONLY FOR 1810 SYRINGES)

Abacus™ shot counter can be purchased pre-installed on Socorex *ultra* 1810 syringe (1810.6.XXXX) or is available as an accessory. It is not compatible with Twin model. It performs precise counts of every injection. A magnet for incrementing manually the counter is available as an accessory.

Counter casing is water-tight to prevent water from entering and damaging the mechanism, for instance in aquaculture applications. **Abacus™ is not autoclavable.**

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Installation (Fig. 18 – 27)

To fit Abacus™ shot counter on Socorex *ultra* 1810 syringe:

1. Unscrew nut and disassemble rod from plunger (Fig. 18).
2. Unscrew socket entirely (Fig. 19).
3. Remove socket from handle (Fig. 20).
4. Remove clip (Fig. 21, ②) in pushed position ①.
5. Assemble spring and socket on magnet-fitted plunger rod supplied (Fig. 22).
6. Insert clip in pushed position (Fig. 23).
7. Place counter holder on socket (Fig. 24, ②), after setting counter orientation among 8 positions available (Fig. 24, ①). Connect to handle ③ (tab fits upper handle recess, Fig. 24, ④).
8. Hold socket while locking ring (Fig. 25).
9. Connect plunger and rod (Fig. 26).
10. Connect body in pushed position. Adjust magnet orientation to that of counter, holding plunger rod motionless (Fig. 27).

Operation (Fig. 28 and 29)

Choose magnet side turning plunger rod:

- Position ①: shot counting.
Position ②: priming, no shot counting.

If position ① was chosen, each plunger activation (shot administered) will increment counter figure.

⚠ **Abacus™ shot counter and magnet-fitted rod are not autoclavable.**

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
WARRANTY

Socorex *ultra* 1810/2810 syringes are guaranteed for a period of one year against any manufacturing or material defect. Cases due to non-respect of manufacturer's instructions are excluded from the warranty. Glass barrel breakage, damages to the plunger or valves due to improper handling or cleaning are not covered by the warranty.

Should regular maintenance not eliminate a detected defect, return the instrument to the dealer from whom it was purchased, after obtaining a return authorisation.

Rinse, dry and decontaminate instrument prior to its return.

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