

# Incubators

PERFECTLY COORDINATED. PERFECTLY CONTROLLED.

INCUBATOR I

CO, INCUBATOR ICOmed

COMPRESSOR-COOLED INCUBATOR ICPeco/ICP

PELTIER-COOLED INCUBATOR IPP

COOLED STORAGE INCUBATOR IPS

100% ATMOSAFE. MADE IN GERMANY.

www.memmert.com | www.atmosafe.net





# Stable. Safe. Sensitive.

Memmert incubators for microbiology. Energy efficient, precise, 100% AtmoSAFE.

Even slight temperature deviations in the working chamber of an incubator may cause a test to fail. For this reason, the heating and control system of Memmert incubators are perfectly adapted to each other. During heating up and cooling down as well as in running operation, all appliances precisely keep the desired parameters within the smallest tolerance limits. Not only at one measuring point, but in the entire working chamber. Each individual Memmert incubator complies with the strict requirements of DIN 12880:2007-05 and is equipped with a maximum of safety functions. Each individual Memmert incubator is 100% AtmoSAFE.



# **INCUBATORS I**

Microbiological tests, colony counts, virology, toxicology

## **CO<sub>2</sub> INCUBATORS ICOmed**

Medical device class IIa for cultivation of cells or tissue, in-vitro fertilisation, gene expression

## COMPRESSOR-COOLED INCUBATORS ICPeco

Microbiological tests, colony counts, virology, toxicology, cultivation above and below room temperature, alternate stability tests

## **COMPRESSOR-COOLED INCUBATORS ICP**

Microbiological tests, colony counts, virology, toxicology, cultivation above and below room temperature, alternate stability tests

## PELTIER-COOLED INCUBATORS IPP

Protein crystallography, microbiological tests, colony counts, virology, toxicology, cultivation above and below room temperature, alternate stability tests

## **COOLED STORAGE INCUBATORS IPS**

Microbiological tests, cultivation above and below room temperature, stability tests

## **ADDITIONAL INFORMATION**

**PAGE 9 - 13** 

**PAGE 4 - 8** 

PAGE 14 - 18

PAGE 19 - 20

PAGE 21 - 25

PAGE 26 - 29

PAGE 30



Incubator IN/INm and IF/IFm with SingleDISPLAY Incubator INplus/INmplus and IFplus/IFmplus with TwinDISPLAY Natural convection or forced air circulation AtmoCONTROL software

Model sizes: 30 / 55 / 75 / 110 / 160 / 260 / 450 / 750 +30 °C to +80 °C

**INCUBATOR I** Memmert incubators I are at home in the world of research, medicine, pharmaceutics and food technology. Organic chamber loads require gentle heating. For this reason, the heating and control system are especially optimised for low temperatures of up to +80 °C. To prevent temperature overshoots, temperature is increased within a very narrow control range and kept exactly at the setpoint value. As required, the models with natural convection or with forced air circulation are available.

4





## As little air circulation as possible in the incubator


Forced air convection may destroy the protective layer from moist air that is generated during incubation over the samples. This would lead to dehydration of the culture. In a Memmert incubator, the perfect combination of all-round surface heating and temperature control system ensures that incubation generally takes place without forced air circulation. Provided the chamber is fully loaded and forced air circulation is required, it can be precisely adjusted in 10 % steps from 0 to 100 %.

## Sterilisation

The chamber of the incubators INplus/IFplus/INmplus and IFmplus, including all installations and sensors, can be sterilised at +160 °C in a 4-hour programme to guarantee optimum hygiene.

## Fresh air is preheated

Temperature deviations caused by fresh air can influence sample characteristics or prolong drying. In Memmert incubators, the fresh air is therefore fed through a pre-heating chamber and seamlessly introduced into the working chamber.



Air supply from outside

#### The incubator Im is a medical device:

The incubator Im is a Class I medical device in accordance with the EU directive 93/42/EEC. In accordance with the intended use incubators INmplus and IFmplus may be used for warming of rinsing solutions and infusions as well as of contrast agents. Incubators INm (with option A6) are intended for heating fango, silicate and APS packs for physical therapy and keeping them warm.

## **INCUBATORS I**


according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

## Standard units are safety-approved and bear the test marks: (EAC not valid for medical devices)

Interior:	Stainless steel, material 1.4301 (ASTM 304) with all-round deep-drawn ribs to integrate the large- area heating with ceramic-metal sheath			
Housing:	Textured stainless steel, rear zinc-plated steel, intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen; inner glass door, outside fully insulated stainless steel door (from size 450 two leaves)			
Fresh air:	Admixture of pre-heated fresh air by electronically adjustable airflap			
Connection:	Mains cable with plug (German type)			
Installation:	4 feet; sizes 450/750 mounted on lockable castors			
Interfaces:	Ethernet USB			

P	
□ LAN □	

USB: only TwinDISPLAY



Model sizes/Desc	cription		30	55	75	110	160	260	450	750	
Stainless steel	Volume	approx. I	32	53	74	108	161	256	449	749	
interior	Width	(A) mm		400		56	50	640	10	40	
	Height	(B) mm	320	400	560	480	720	800	720	1200	
	Depth (less 39 mm for fan)	(C) mm	250	33	30	40	00	500	6	00	
	Max. number of grids/shelves	number	3	4	6	5	8	9	8	14	
	Max. loading per grid/shelf	kg			2	20			3	0	
	Max. loading of chamber	kg	60	80	120	175	210		300		
	Max. loading per slide-in drip tray	kg		1,5		3	3	4	1	3	
	Max. loading per bottom drip tray	kg		1,5			3	4		3	
Textured stainless	Width	(D) mm		585		74	45	824	12	24	
steel exterior	Height (size 450, 750 with castors)	(E) mm	704	784	944	864	1104	1183	1247	1726	
	Depth (without door handle), door handle + 56 mm	(F) mm	434	5	14	58	34	684	78	34	
Standard	Stainless steel grids, electropolished	number		1				2			
equipment	Standard works calibration certificate (measuring point chamber center)	°C	+37								
Temperature	Working temperature range	°C	min. 5 (IN/INplus/INm/INmplus) 10 (IF/IFplus/IFm/IFmplus) above ambient temperature up to +80								
	Setting temperature range	°C	+20 to +80								
	Setting accuracy	°C	0.1								
Further data	Electrical load at 230 V, 50/60 Hz	approx. W	1600	1000	1250	1400	1600	1700	1800	2000	
	Electrical load at 115 V, 50/60 Hz	approx. W	800	900					1500	1800	
Packing data	Net weight	approx. kg	48	57	66	76	96	110	161	217	
J	Gross weight (packed in carton)	approx. kg	64	76	85	101	122	161	227	288	
	Width	approx. mm	660	7	30	83	30	930		30	
	Height	approx. mm	890	950	1130	1050	1300	1380	1440	1910	
	Depth	approx. mm	650	6	70	80	00	930	10	50	
Order No. Incubators		IN30	IN55	IN75	IN110	IN160	IN260	IN450	IN750		
I = Incubator		IN30plus	IN55plus	IN75plus	IN110plus	IN160plus	IN260plus	IN450plus	IN750plus		
N = Natural cor	nvection		IF30	IF55	IF75	IF110	IF160	IF260	IF450	IF750	
F = Forced convection m = Medical device			IF30plus	IF55plus	IF75plus	IF110plus	IF160plus	IF260plus	IF450plus	IF750plus	

plus = Model with TwinDISPLAY

Options	30	55	75	110	160	260	450	750	
Voltage 115 V, 50/60 Hz				X	2				
Extended overtemperature protection by additionally integrated Pt100 sensor for independent temperature monitoring for models with SingleDISPLAY	A6								
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids		-						(1	
Interior lighting for observing the load				F	10				
Interior socket (can only be ordered with limited temperature range - max. +70 °C) ampacity 230 V, 2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (only with SingleDISPLAY), (option A8 necessary)				R	3				
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap and silicone stopper, standard positions	F2								
Entry port, 23 mm clear diameter, left					4				
can be closed by flap, in special right positions (please state location) rear				F	5 6				
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)				۵	06				
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location) Entry port, 57 mm clear diameter, can be closed by flap, in special positions					7				
in the back wall (please, state location) Entry port, 100 mm clear diameter, can be closed by flap, in special				F	8				
positions in the back wall (please, state location)	-	-				-9			
4 - 20 mA current loop interface (0 to +90 °C = 4 - 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY) – price per sensor					/3 /6				
Fan speed monitoring with switching-off the heating and with alarm in case of failure - only for IFplus/IFmplus	V4								
Works calibration certificate for 3 temperatures: +37 °C, +52 °C, +70 °C	D00126								
Works calibration certificate for one (freely selectable) temperature value according to customer specification	D00109								
Door with lock and key (safety lock)					6				
Door hinged on the left Potential-free contact (24 V/2 A) with socket, according to NAMUR NE 28			В					-	
for external monitoring (indicates when setpoint is reached) Potential-free contact for combination error message (e.g. supply failure,	Н5 Н6								
sensor fault, fuse)									
Potential-free contact (24 V/2 A) 2 contacts with socket to NAMUR NE 28, for signal generation, controlled by programme segment, for free- selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY	ts H72								
Process-dependent programmable door lock (only for units with TwinDISPLAY)	D4								
Door-open-recognition (only for units with TwinDISPLAY)	V5								
Flexible Pt100 for positioning in chamber or in load with socket, 4-pin, according to NAMUR NE 28, for external temperature recording (load temperature) max. 3 sensors	H4								
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software	e H8								
MobileALERT, notification by SMS in case of any error or alarm of the device. Requires option H6				C	3				
Temperature restriction; Temperatures: +50, +55, +60, +65, +70 or +75 °C (Please, indicate upon ordering)	A8								
Fresh-air filter (filtration efficiency 80 %) mounted at the appliance bottom (for IF/IFplus/IFm/IFmplus). For sizes 30 – 260 castor frame R9 or subframe necessary				F	8				
Castor frame (2-part), height 140 mm			R	9				-	

Accessories	30	55	75	110	160	260	450	750
Stainless steel grid, electropolished	E28884	E28884 E20164 E20165 E288			E28891	E20	E20182	
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		- E29767 E29766			E29766	B32	190	
Perforated stainless steel shelf	B29727	B03	916	BOC	325	B29725	B00	328
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber				-			B32191	
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	E02070	E02	072	E02	073	E29726	E02	075
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)				-			B32	763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	B04356	B04	358	B04	359	B29722	B04	362
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)				-			B34	055
Wall bracket for wall mounting	B29755	B29756	B29757	B29758	B29759		-	
Guarantee extension by 1 year			GA1Q5				GA2Q5	
USB-Ethernet adapter				E06	5192			
Ethernet connection cable 5 m for computer interface				E06	5189			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID- programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)				B33	3170			
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number	B33172							
Set of height adjustable feet (4 pcs)			B29	768				-
Stacking set (4 pcs) for stacking of appliances of same size		B29	744				-	
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)				B29	9718			
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)				B29	9719			
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	B29738	B29740	B29742
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29733	B29735	B29737	B29739	B29741	B29743
Subframe, adjustable in height (size 30 to 75: height 600 mm, size 110 to 450: height 500 mm)	B29745	B29	747	B29	749	B29751	B29753	-
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746	B29	748	B29	750		-	
Subframe, adjustable in height, height 130 mm, for example for units with fresh air filter	B33657	B33	659	B33	661	B33664		-
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1							
Integration of one additional unit (up to max. 15 units) into an already existent FDA- software licence (only for units with TwinDISPLAY)	FDAQ2							
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer	D00124							
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 - 1060) to DIN 12880:2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer site on demand (GER, AT, CH only)	D00125 D00127							



CO<sub>2</sub> Incubator ICOmed with TwinDISPLAY Software AtmoCONTROL

Model sizes: 50 / 105 / 150 / 240+18 °C to +50 °C Humidity 40 to 97 % rh CO<sub>2</sub> concentration 0 to 20 % O<sub>2</sub> concentration 1 to 20 %

**CO**<sub>2</sub> **INCUBATOR ICOmed** Safety at all times. When it comes to safety and user friendliness, the highly modern  $CO_2$  incubator ICOmed is the perfect solution: Thanks to the battery-buffered ControlCOCKPIT, the operating display, logging and  $CO_2$  control remain fully functional even when there is a power failure. All parameters are logged in accordance with the FDA and, when individually adjusted ranges for  $CO_2$ ,  $O_2$ , temperature and humidity are exceeded, notifications can be sent to a mobile phone in addition to an alarm.

The control technology is so finely tuned that the setpoint temperature is guaranteed to be reached without temperature overshoots. With its rounded corners, the interior is easy to clean and can be sterilised for 60 minutes at 180 °C (including all sensors).

All ICOmed models are classified as medical product class IIa.





# Unrivalled user friendliness

All parameters can be set easily and intuitively both with the ControlCOCKPIT or the AtmoCONTROL software. The shutter box can be opened, allowing fast access to controls. Maintenance is possible even if the appliances are stacked. The appliance has USB and Ethernet connections as well as a data logger with a ten-year storage capacity. Data can be read and programmes can be transferred by remote access.

# Minimising vaporisation and condensation

The active humidity control minimises vaporisation in the interior and ensures short recovery times after the door has been opened. Together with the heating of the interior from all six sides including the heated inner glass door, it prevents the dangerous formation of condensation and offers maximum protection for cell and tissue cultures. The turbulence-free chamber ventilation ensures a constant and uniform atmosphere.



# IVF module for models ICO50med/ICO105med

In order to keep vaporisation, condensation and recovery times at a minimum during in vitro fertilisation, the Petri dishes are cultivated in separate slide-in units. The slide-in units in the optional IVF module can be pulled out easily and with low vibration and are equipped with a pull-out lock.



### The CO<sub>2</sub> Incubator ICOmed is a medical device:

Memmert subjected its  $CO_2$  incubator ICOmed to a comprehensive medical device evaluation. Every Memmert  $CO_2$  incubator ICOmed is classified as a Class IIa medical device. The ICOmed is intended for the creation and maintenance of constant environmental conditions for application in the field of in vitro fertilisation (IVF), especially for the incubation of oocytes, spermatozoa and zygotes in special culture dishes for IVF application as well as for gene expression and the biosynthesis of RNA and proteins. The CE label on the appliances includes the mark 0197, denoting TRLP – TÜV Rheinland as the notified body.

56

## **CO<sub>2</sub> INCUBATORS ICOmed**

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

Standard units are safety-approved and bear the test marks:  $\bigcup_{U \in U} \bigcup_{U \in U} \bigcup_{U$ 

 Interior:
 Stainless steel, material 1.4301 (ASTM 304), deep-drawn, seamlessly welded


 Housing:
 Textured stainless steel, rear zinc-plated steel, intuitively operated TwinDISPLAY (TFT colour display) with touchsceen; fully insulated stainless steel door and heated inner glass door

 Automatic
 Humidity and CO<sub>2</sub> sensor sterilised inside the CO<sub>2</sub> incubator

 Connection:
 Mains cable with plug (German type)

 Installation:
 4 adjustable feet

 Interfaces:
 Ethernet



Model sizes/Des	ription		50	105	150	240	
Stainless steel	Volume	approx. I	56	107	156	241	
interior	Width	(A) mm	400	0 560		600	
	Height	(B) mm	425	480	700	810	
	Depth (less 35 mm for fan)	(C) mm	330	40	00	500	
	Max. number of perforated shelves	number	5	6	10	12	
	Max. loading per perforated shelf	kg		1	15		
	Max. loading of chamber	kg	75	90	120	140	
Textured stainless	Width	(D) mm	559	7	19	759	
steel exterior	Height (variable through adjustable feet)	(E) mm	795	850	1070	1180	
	Depth (without door handle, depth of door handle 56 mm)	(F) mm	521		91	691	
	Fully insulated heated stainless steel door	()					
	Additional heated inner glass door						
Ctore dowd			1		2		
Standard equipment	Stainless steel shelves, perforated	number	I		2		
- 4	Stainless steel water dish (not applicable with option K7) Works calibration certificate (measuring point chamber centre) at +37 °C, 5 % CO <sub>2</sub> for standard units	number			•		
	Works calibration certificate at 37 °C, 5 % CO <sub>2</sub> , 90 % rh and 10 % O <sub>2</sub> (requires option K7 and option T6); standard equipment for units with O <sub>2</sub> control			(	Đ		
	Works calibration certificate at 37 °C, 5 % $CO_2$ and 90 % rh (requires option K7); standard equipment for units with active humidity control		•				
	CO <sub>2</sub> connection set: hose with coupling and clamp				•		
	Standard sterilisation programme (without removing the sensors), humidity and $CO_2$ sensor sterilised inside the $CO_2$ incubator		60 minutes at 180 °C				
	Membrane filter (in order to remove impurities and pollutants, all incoming gases pass through a membrane filter before they reach the chamber)				•		
Temperature	Working temperature range	°C	5 above	e ambient tei	mperature u	o to +50	
	Setting temperature range	°C	+18 to +50				
	Setting accuracy	°C	0.1				
	Temperature fluctuations with time at 37 °C (to DIN 12880:2007-05)	К	+/- 0.1				
	Temperature variation in chamber at + 37 °C (to DIN 12880:2007-05)	К		+/-	- 0.3		
Humidity	Humidity limitation thanks to a Peltier element; when water dish is full and inserted, the Peltier element limits the value of relative humidity in the interior to 93 % rh +/- 2.5 %				D		
	Setting range active humidity control (with option K7)	% rh	40 to 97 and rh-Off				
	Setting accuracy	% rh		0	.5		
CO <sub>2</sub> / O <sub>2</sub>	Digital electronic CO <sub>2</sub> control with dual beam NDIR system, with auto-diagnostic system and acoustic fault indication, barometric pressure compensation				•		
	Setting range CO <sub>2</sub>	% CO <sub>2</sub>		0 to	o 20		
	Variation in time CO <sub>2</sub>	% CO <sub>2</sub>	+/- 0.2				
	Setting accuracy CO <sub>2</sub>	% CO <sub>2</sub>	0.1				
	Setting range O <sub>2</sub>	% O <sub>2</sub>			o 20		
	Setting accuracy O <sub>2</sub>	% 0 <sub>2</sub>			0.1		
Further data	Electrical load at 230/115 V, 50/60 Hz	approx. W	1000	1500		000	
Packing data	Net weight	approx. kg	55	75	90	110	

Model sizes/Desc	ription			50	105	150	240	
Packing data	Gross weight (packed in c	arton)	approx. kg		100	116	145	
	Width		approx. mm		80		840	
	Height Depth		approx. mm approx. mm		1030 80	1250	1360 900	
Order No. CO <sub>2</sub> Incu			αρριολ. ΠΠ	ICO50	ICO105	ICO150	IC0240	
Options			50	105	150		240	
/oltage 115 V, 50/60	) Hz		50	X	_		240	
Battery-buffered Con	trolCOCKPIT: Uninterrupted ocumentation of all parame	supply for the entire display unit (ControlCOCKPIT) and eters even when there is a power failure. The $\rm CO_2$		C				
wo gas connections	, ,	ors for automatic switch-over of gas cylinders; incl. two $\mathrm{CO}_2$		T	1			
lectropolished interi	or			T	2			
auto-diagnostic syste while avoiding conde	em ensures even more rapio ensate formation. Humidity	d dehumidifying (40 - 97 % rh), incl. digital indication and I reaching of set humidity and very short recovery times supply with distilled water (from an external tank) by a nerating hotsteam, dehumidifying via sterile filter		K	7			
Control of oxygen cor requires option K7).	ncentration by $N_2$ inlet; adjuined in $N_2$ connection set: hose the set in $N_2$ connection set in $N_2$ connection set in $N_2$ connection set in $N_2$ connection set in $N_2$ set in	istment range 1 % up to 20 % ${\rm O_2}$ ; setting accuracy 0.1 % e with coupling and clamp		T	5			
Peltier cooling unit er 85 °C	nables a working temperati	ure of 37 °C even at higher ambient temperatures of up to		K	5			
apacitive humidity	sensor for measuring and d	isplaying the relative humidity		K	6			
ntry port (silicone), 4 entre right; not avai option K6)	40 mm clear diameter, mois lable for ICO50med with ac	ture tight, can be closed by silicone stopper, at the back, tive humidity control (option K7) or humidity display	F7					
nner door with partit as 2/3/4 partitioned		be used in connection with option B8); size 105/150/240	- K4					
- 20 mA current loc	p interface	Temperature controller, actual value (0 to +70 °C = $4 - 20$ mA)		V.	3			
		Humidity controller, actual value (0 to 100 % $rh = 4$ - 20 mA) (requires option K7 or K6)		V	7			
		$CO_2$ controller, actual value (0 to 25 % $CO2 = 4 - 20$ mA)	V9					
		$O_2$ controller, actual value (0 to 25 % $O2 = 4 - 20$ mA) (requires option T6)	V1					
Vorks calibration cer ertificates upon requ		$\%$ CO $_{\rm 2}$ (measured at +37 °C) special works calibration	D00106					
Vorks calibration cer ustomer specificatio	tificate for one (freely selec n (requires K7)	table) temperature, humidity and $CO_2$ value according to	D00131					
Vorks calibration cer o customer specifica		table) temperature, humidity, $CO_2$ and $O_2$ value according	D00143					
tart-up of ICOmed ir liscount	ncubators and brief training	(D, A, CH only), through Memmert service, not subject to	К9					
oor hinged on the le			B8					
Potential-free contact (24 V/2 A) with socket to NAMUR NE 28 for external monitoring (indicates when set points of temperature and $CO_2$ are reached)			H5					
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)			H6					
MobileALERT, notification by SMS in case of any error or alarm of the device. Requires option H6 MobileALERT for up to 4 alarm notifications; temperature and CO <sub>2</sub> alarm (standard), additionally humidity			C3					
viobileALERT for up t alarm (when equippe	to 4 alarm notifications; ten ed with option K7) and $O_2$ a	perature and $CO_2$ alarm (standard), additionally humidity larm (when equipped with option T6)		C	4			
Accessories				50	105	150	240	
erforated stainless s	teel shelf			E3516	0 E3	7418	E3515	
Vater dish					BB	8737		
ubframe (622 mm h	nigh) adjustable in height (s	izes 150/240: should not be used for 2 stacked units)		B3350	4 B3	3505	B3350	
ubframe (130 mm h		n combination with the corresponding stacking sets for stacked	d appliances	B3350	7 B3	3508	B3350	

HEPA-filter for chamber (filter class E11) according to EN 1822, packed in sterile condition, incl. fixing unit

Central water supply, with filter cartridges for connection to the domestic water supply, only in combination with option K7. Product

Central water supply, without filter cartridges for connection to the domestic water supply (only for demineralised water in accordance with VDE 0510/DIN EN 50272), only in combination with option K7. Product information on demand

IVF-module for ICO50med: patented consisting of 6 slide-in units, a total of 12 special racks with indentations for 12 Petri dishes (60 mm diam.) resp. 24 Petri dishes (35 mm diam.), 2 racks with indentations for 3 special media tubes each; racks with indentations for 4-well dishes on demand; only for ICO50med with the options K7 and F7; works calibration certificate (measuring point chamber centre) at +37 °C, 5 %, 6 % and 7 % CO<sub>2</sub> as well as 90 % rh; 5 % O<sub>2</sub> for IVF unit equipped with option T6

N<sub>2</sub> pressure reducing valve to DIN EN ISO 2503, incl. gas cylinder monitor (requires option T6)

CO2 pressure reducing valve to DIN 8546, incl. gas cylinder monitor

Celltron benchtop shaker (not subject to discount) - accessories upon request

B38739

E02087

E06162

ZWVR6

ZWVR7 GA3Q5

E06724

B44128

information on demand

Guarantee extension by 1 year

Accessories	50	105	150	240
IVF-module for ICO105med: patented, consisting of 8 slide-in units, a total of 16 special racks with indentations for 16 Petri dishes (60 mm) resp. 32 Petri dishes (35 mm diam.), 2 racks with indentations for 3 special media tubes each; racks with indentations for 4-well dishes on demand; only for ICO105med with the options K7 and F7; works calibration certificate (measuring point chamber centre) at 37 °C, 5 %, 6 % and 7 % CO <sub>2</sub> as well as 90 % rh; 5 % O <sub>2</sub> for IVF unit equipped with option T6	-	B42398		-
Holder for Petri dishes round (only in combination with IVF-module)	E37	026		-
Holder for Petri dishes square (only in combination with IVF-module)	E37	308		-
Holder for test tube (only in combination with IVF-module)	E37	069		-
Magnetic foil, can be labelled with a non-permanent board marker (only in combination with IVF-module)	E36	651		-
USB-Ethernet adapter		E06	192	
Ethernet connection cable 5 m for computer interface		E06	189	
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170			
Stacking set (4 pcs) for stacking of appliances of same size	B29	9744		-
Stacking set (consisting of stacking corners, one connecting plate for the rear, two wall brackets) for stacking of two units of same size		-	B42114	B42115
FDA confroming software AtmoCONTROL (FDA edition). Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit. Respective IQ/OQ documents available in German and English language (without surcharge)	tmoCONTROL (FDA edition). Meets the requirements for the use of electronically stored data sets and down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the ye IO/OO documents available in German and English language (without surcharge)		AQ1	
Integration of additional units (up to max. 15 units) into an already existent FDA-software licence	FDAQ2			
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer	D00124			
IQ/OQ document with device-specific works test data for one free-selectable $CO_2$ , humidity and temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05, PQ check list as support for validation by customer (a free-selectable humidity value is only possible with option K7). Price for validation at customer site on demand (GER, AT, CH only)	D38897			
IQ/OQ document with device-specific works test data for one free-selectable $CO_2$ and temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05, PQ check list as support for validation by customer. Price for validation at customer site on demand (GER, AT, CH only)	D38898			
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand B04714			714	



CO<sub>2</sub>-cooled incubator ICPeco with TwinDISPLAY AtmoCONTROL software

Model sizes: 110 / 260 / 450 / 750 (ICP110eco available from Q2 2019) -12 °C to +60 °C

**COMPRESSOR-COOLED INCUBATOR ICPeco** These environmentally-friendly cooled incubators are cooled with climate-neutral CO<sub>2</sub>. Thanks to this refrigerant's excellent thermodynamic properties and the finely adjusted control technology, an ICPeco is both powerful and high-precision. Without critical temperature overshoots, it keeps the temperatures exactly at the setpoint.






## Refrigerant CO<sub>2</sub> is climate-neutral

A CO<sub>2</sub>-cooled incubator ICPeco is in every respect positive for the ecological balance of a laboratory. Legal restrictions for use are completely excluded in the future, as the refrigerant  $CO_2$  (R744), unlike fluorine-based refrigerants, has no greenhouse gas reduction potential. It is a by-product of industrial processes, which is why far less energy is used for its manufacture than for synthetic, fluorinated refrigerants. R744 is neither flammable nor toxic, does not cause ozone depletion in the atmosphere and does not require disposal or recycling.

# Refrigerant CO<sub>2</sub> ensures better cooling performance

The contribution to process optimisation is also impressive. An ICPeco is virtually maintenancefree and extremely powerful. Compared to appliances with R134a as refrigerant, it has faster temperature change rates during cooling-down. Memmert cooled incubators ICP with refrigerant R134a will be available in parallel for a transitional period.




On average 20 % faster temperature change rates with a CO<sub>2</sub>-cooled compressor (measurement ICP260eco at ambient temperature +22 °C according to IEC 60068-3-5)





Cooling and heating units are situated outside the working chamber inside the air jacket temperature control system surrounding the entire chamber interior ensuring quick and precise temperature control. The motor-driven forced air circulation, adjustable in 10 % steps via the ControlCOCKPIT ensures optimum temperature distribution.





ICP air jacket temperature control system

# Integrated energy saving function


The cooling unit works extremely energy-efficient because the heating is completely switched off in cooling mode. An intelligent DEFROST function enables defrosting as required.


## COMPRESSOR-COOLED INCUBATORS ICPeco

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

# Standard units are safety-approved and bear the test marks: $\mathsf{C} \in \mathsf{E} \mathsf{F} \mathsf{E} \mathsf{F} \mathsf{E}$

Interior:	Stainless steel, material 1.4301 (ASTM 304)				
Housing:	Textured stainless steel, rear zinc-plated steel, intuitively operated TwinDISPLAY (TFT colour display) with touchscreen; inside glass door, outside fully insulated stainless steel door (from size 450 two leaves)				
Connection:	Mains cable with plug (German type)				
Installation:	Mounted on lockable castors				
Interfaces:	Ethernet USB				





Model sizes/Desc	ription		110	260	450	750
Stainless steel	Volume	approx. I	108	256	449	749
interior	Width	(A) mm	560	640	10	)40
	Height	(B) mm	480	800	720	1200
	Depth (less 33 mm for fan)	(C) mm	400	500	6	00
	Max. number of grids/shelves	number	5	9	8	14
	Max. loading per grid/shelf	kg	2	0	3	80
	Max. loading of chamber	kg	150		200	
	Max. loading per slide-in drip tray	kg	3	4		8
	Max. loading per bottom drip tray	kg	3	4		8
Textured stainless	Width	(D) mm	745	824	12	24
steel exterior	Height (with castors)	(E) mm	1233	1552	1613	1950
	Depth (without door handle), door handle + 56 mm	(F) mm	584	684	7	84
Standard	Stainless steel grids, electropolished	number	mber 2			
equipment	Standard works calibration certificate (measuring point chamber center)	°C		+10 a	nd +37	
Temperature	Working temperature range (not suitable for long-term storing at sub-zero temperatures. During permanent operation, the glass door may ice over)	°C		-12 t	o +60	
	Setting temperature range	°C		-12 t	io +60	
	Setting accuracy	°C		C	).1	
Further data	Electrical load at 230 V, 50 Hz	approx. W		12	200	
Packing data	Net weight	approx. kg	118	162	222	254
	Gross weight (packed in carton)	approx. kg	146	219	287	324
	Width	approx. mm	880	930	13	30
	Height	approx. mm	1410	1760	1700	2150
	Depth	approx. mm	810	930	10	)50
Order No. Compre	essor-Cooled Incubators					

(ICP110eco available from Q2 2019)

ICP110eco ICP260eco ICP450eco ICP750eco

Options		110	260	450		750
Chamber modification for the application of reinfo grids (bearing rails mounted in the working cham reinforced grids	-			K1		
Interior socket, ampacity 230 V/2.2 A, can be swit individually, moisture tight IP68	ched off with the On/Off switch, cannot be switched		R3			
Entry port, 23 mm clear diameter, for introducing	left centre/centre		FO			
connections at the side, can be closed by flap and silicone stopper, standard positions	left centre/top		F1			
	right centre/top	-		F3		
Entry port (silicone), 40 mm clear diameter, moistı (please, state location)	ure tight, can be closed by silicone stopper, at the back		F7			
4 - 20 mA current loop interface	Temperature controller, actual value (-20 to +70 °C = 4 - 20 mA) Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max.		V3 V6			
	3) - price per sensor (-20 to +70 °C = 4 - 20 mA)					
Fan speed monitoring with switching off the heat	ing and with alarm in case of failure		V4			
Works calibration certificate for 3 temperatures: 0			D0013	0		
Works calibration certificate for one (freely selecta specification	ble) temperature value according to customer		D0010	9		
Door with lock and key (safety lock)			B6			
Door hinged on the left		B8			-	
Potential-free contact (24 V/2 A) with socket, acco when setpoint is reached)	ording to NAMUR NE 28 for external monitoring (indicates		H5			
Potential-free contact for combination error messa	ge (e.g. supply failure, sensor fault, fuse)		H6			
Potential-free contact (24 V/2 A) with socket to NAMUR NE 28, for signal generation, controlled b programme segment, for free-selectable functions to be activated (e.g. activation of audible and visu signals, exhaust motors, fans, stirrers, etc.)	5	Н72				
Process-dependent programmable door lock			D4			
Door-open-recognition		V5				
Flexible Pt100 for positioning in chamber or in loa external temperature recording (load temperature	d with socket, 4-pin, according to NAMUR NE 28, for e) max. 3 sensors		H4			
Flexible Pt100 temperature sensor, positioned flex measurement (up to 2 additional sensors are poss indicated on the display, recorded in the integral of AtmoCONTROL software	sible). The measured temperature can, if required, be	H8				
MobileALERT, notification by SMS in case of any en	rror or alarm of the device. Requires option H6		G			
Accessories			110	260	450	750
Stainless steel grid, electropolished			E20165	E28891	E20	182
	olished, max. loading 60 kg; from size 450 with guide bars ar loading of chamber	nd fixing screws (only ir		E29766		190
Perforated stainless steel shelf					BOC	)328
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber			B00325	B29725		191
tainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1			E02073	E29726	E02	075
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)				-		2763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1						
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)				-		1055
USB-Ethernet adapter				E06	192	
Ethernet connection cable 5 m for computer interface				E06182		
	ed authorisation licence (User-ID-programme) on Memory-st	tick, prevents undesired	I		170	
Flush-fit unit (stainless steel frame covering gap b			B29734	B29738	B29740	B29742
Flush-fit unit (stainless steel frame covering gap b	1 5.0		B29735	B29739	B29741	B2974


Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots B29735 B29739 B29741 B29743 FDA confroming software AtmoCONTROL (FDA edition). Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit. Respective IQ/OQ documents available in German and English language (without surcharge) FDAQ1 Integration of additional units (up to max. 15 units) into an already existent FDA-software licence FDAQ2 IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer D00124 IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer site on demand (GER, AT, CH only) D00127


## **COMPRESSOR-COOLED INCUBATORS ICP**

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

# Standard units are safety-approved and bear the test marks: $\mathsf{C} \in \mathsf{E}\mathsf{R}\mathsf{I}$

Interior:	Stainless steel, material 1.4301 (ASTM 304)				
Housing:	Textured stainless steel, rear zinc-plated steel, intuitively operated TwinDISPLAY (TFT colour display) with touchscreen; inside glass door, outside fully insulated stainless steel door (from size 450 two leaves)				
Connection:	Mains cable with plug (German type)				
Installation:	Mounted on lockable castors				
Interfaces:	Ethernet USB				





Model sizes/Deso	ription		110	260	450	750
Stainless steel	Volume	approx. I	108	256	449	749
interior	Width	(A) mm	560	640	10	40
	Height	(B) mm	480	800	720	1200
	Depth (less 33 mm for fan)	(C) mm	400	500	6	00
	Max. number of grids/shelves	number	5	9	8	14
	Max. loading per grid/shelf	kg	2	0	3	0
	Max. loading of chamber	kg	150		200	
	Max. loading per slide-in drip tray	kg	3	4		8
	Max. loading per bottom drip tray	kg	3	4		8
Textured stainless	Width	(D) mm	745	824	12	24
steel exterior	Height (with castors)	(E) mm	1233	1552	1613	1950
	Depth (without door handle), door handle + 56 mm	(F) mm	584	684	7	84
Standard	Stainless steel grids, electropolished	number		-	2	
equipment	Standard works calibration certificate (measuring point chamber center)	°C		+10 ar	nd +37	
Temperature	Working temperature range (not suitable for long-term storing at sub-zero temperatures. During permanent operation, the glass door may ice over)	°C		-12 to	o +60	
	Setting temperature range	°C		-12 to	o +60	
	Setting accuracy	°C		0	.1	
Further data	Electrical load at 230/115 V, 50/60 Hz	approx. W		12	00	
Packing data	Net weight	approx. kg	113	157	217	249
-	Gross weight (packed in carton)	approx. kg	141	214	282	319
	Width	approx. mm	880	930	13	30
	Height	approx. mm	1410	1760	1700	2150
	Depth	approx. mm	810	930	10	50
Order No. Compre	essor-Cooled Incubators		ICP110	ICP260	ICP450	ICP750

Options		110	260	450		750
Voltage 115 V, 50/60 Hz			X2			
Chamber modification for the application of reinfor grids (bearing rails mounted in the working chamb reinforced grids	ced perforated stainless steel shelves or stainless steel per) - includes replacement of standard grids by	- K1				
Interior socket, ampacity 230 V/2.2 A, can be swite individually, moisture tight IP68	hed off with the On/Off switch, cannot be switched		R3			
Entry port, 23 mm clear diameter, for introducing	left centre/centre		FO			
connections at the side, can be closed by flap and silicone stopper, standard positions	left centre/top		F1			
	right centre/top	- F3				
Entry port (silicone), 40 mm clear diameter, moistu (please, state location)	re tight, can be closed by silicone stopper, at the back		F7			
4 - 20 mA current loop interface	4 - 20 mA current loop interface (-20 to +70 °C = 4 to 20 mA)		V3			
	Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 3) - price per sensor (-20 to +70 °C = 4 - 20 mA)		V6			
Fan speed monitoring with switching off the heating	ng and with alarm in case of failure		V4			
Works calibration certificate for 3 temperatures: 0 °	°C, +37 °C, +60 °C		D00130	)		
Works calibration certificate for one (freely selectable specification	le) temperature value according to customer	D00109				
Door with lock and key (safety lock)		B6				
Door hinged on the left		B8 -				
Potential-free contact (24 V/2 A) with socket, account when setpoint is reached)	rding to NAMUR NE 28 for external monitoring (indicates		H5			
Potential-free contact for combination error message	ge (e.g. supply failure, sensor fault, fuse)		H6			
Potential-free contact (24 V/2 A) with socket to NAMUR NE 28, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visua signals, exhaust motors, fans, stirrers, etc.)	2 contacts		H72			
Process-dependent programmable door lock			D4			
Door-open-recognition			V5			
· · ·	l with socket, 4-pin, according to NAMUR NE 28, for max. 3 sensors		H4			
Flexible Pt100 temperature sensor, positioned flexi measurement (up to 3 additional sensors are poss indicated on the display, recorded in the integral d AtmoCONTROL software	ble). The measured temperature can, if required, be	H8				
MobileALERT, notification by SMS in case of any en	or or alarm of the device. Requires option H6		G			
Accessories			110	260	450	750
Stainless steel grid, electropolished			E20165	E28891	E20	182
Additional reinforced stainless steel grid, electropo connection with option K1). Please consider max. I	lished, max. loading 60 kg; from size 450 with guide bars ar oading of chamber	nd fixing screws (only i	n E29767	E29766	B32	190
Perforated stainless steel shelf		B00325	B29725	B00	328	

Perforated stainless steel shelf	B00325	B29725	BOC	)328
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber		-	B32	2191
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	E02073	E29726	E02	075
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)		-	B32	2763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	B04359	B29722	B04	362
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)	- B34055			1055
USB-Ethernet adapter		E06	192	
Ethernet connection cable 5 m for computer interface	E06189			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number	B33170			
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29734 B29738 B29740 B2		B29742	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29735	B29739	B29741	B29743
FDA confroming software AtmoCONTROL (FDA edition). Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit. Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1			
Integration of additional units (up to max. 15 units) into an already existent FDA-software licence	FDAQ2			
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer		D00	124	
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880:2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer site on demand (GER, AT, CH only)	t D00127			



Peltier-cooled incubator IPP with SingleDISPLAY Peltier-cooled incubator IPPplus with TwinDISPLAY AtmoCONTROL software

Model sizes: 30 / 55 / 110 / 260 / 400 / 750 / 1060 0 °C to +70 °C

**PELTIER-COOLED INCUBATOR IPP** Heating and cooling seamlessly with one system thanks to Peltier technology. In this respect, cooled incubators IPP not only contribute to climate protection, but it also achieves an additional decrease in operating costs of up 90 % compared to compressor technology. This perfect development from the environmentally friendly and energy-saving heating/cooling technology by Memmert convinces by outstanding control precision and extremely small fluctuations.





#### Extremely quiet and vibration-free


The fact that no compressor is required saves space and brings peace and quiet to the laboratory. As Peltier-cooled incubators IPP are almost vibration-free, they can also be applied in entomology. If defined humidity is also required, an alternative would be the constant climate chamber HPP, which is also equipped with Peltier technology.

# No condensation in the interior chamber

Due to the closed Peltier cooling system, no outside air is exchanged. Physically derived, unavoidable formation of condensation during the cooling process does not take place in the interior chamber but on the outside heat sink. In addition, the in the Peltier elements integrated fans ensure a rapid transport of energy as well as an optimal temperature distribution.

# Energy-saving heating/cooling technology combination

In contrast to compressor systems, Peltier technology is particularly economical at temperatures close to the ambient temperature, since energy is only required during heating or cooling. Therefore heating and cooling function are particularly precisely adjusted to each other.






## PELTIER-COOLED INCUBATORS IPP

## according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

# Standard units are safety-approved and bear the test marks: $\mathsf{C} \in \mathsf{E} \mathsf{F} \mathsf{E} \mathsf{F} \mathsf{E}$

Interior:	Stainless steel, material 1.4301 (ASTM 304), deep- drawn	(93) D (93)	3)
Housing:	Textured stainless steel, rear zinc-plated steel, intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen		
Double doors:	Outside stainless steel, fully insulated, inside glass (size 750 and 1060 two leaves)		319
Connection:	Mains cable with plug (German type)		
Installation:	4 feet; sizes 400 to 1060 mounted on lockable castors		
Interfaces:	Ethernet LAN USB USB: only for TwinDISPLAY		30



Model sizes/Desc	ription		30	55	110	260	400	750	1060	
Stainless steel	Volume	approx. I	32	53	108	256	384	749	1060	
interior	Width	(A) mm	4	00	560	64	40	10	1040	
	Height	(B) mm	320	400	480	800		1200		
	Depth (less 10 mm for fan – Peltier)	(C) mm	250	330	400	50	00	600	850	
	Max. number of grids/shelves	number	3	4	5	9		14		
	Max. loading per grid/shelf	kg			20			30	20	
	Max. loading of chamber	kg	60	80	150		2	00		
	Max. loading per slide-in drip tray	kg	1	,5	3		4		3	
	Max. loading per bottom drip tray	kg	1	,5	3		4		3	
Textured stainless	Width	(D) mm	5	85	745	8	24	1224		
steel exterior	Height (sizes 400, 750, 1060 with castors)	(E) mm	704	784	864	1183	1720	17	26	
	Depth (without door handle), door handle + 56 mm	(F) mm	506	586	656	7	56	856	1107	
Standard	Stainless steel grids, electropolished	number		1 2						
equipment	Standard works calibration certificate (measuring point chamber center)	°C				+10 and +3	7			
Temperature	Working temperature range without light	°C		0 (at least 20 below ambient temperature) to +70						
	Working temperature range with light	°C		- +10 to +40				-		
	Setting temperature range	°C		0 to +70						
	Setting accuracy	°C				0.1				
Further data	Electrical load at 230/115 V, 50/60 Hz	approx. W	140	275	550	820	1100	1300	1500	
	Peltier elements in the rear	number		1	2	3	5		5	
Packing data	Net weight	approx. kg	40	52	78	114	157	230	255	
-	Gross weight (packed in carton)	approx. kg	56	71	103	165	210	301	419	
	Width	approx. mm	660	730	830	93	30	1330	1370	
	Height	approx. mm	890	950	1050	1380	1930	1910	1970	
Depth		approx. mm	650	670	800	93	30	1050	1300	
Order No. Peltier-	Cooled Incubators		IPP30	IPP55	IPP110	IPP260	IPP400	IPP750	IPP1060	
IPP = Peltier-Cool			IPP30plus	IPP55plus	IPP110plus	IPP260plus	IPP400plus	IPP750plus	IPP1060plus	

plus = Model with TwinDISPLAY

Options	30	55	110	260	400	750	1060
Voltage 115 V, 50/60 Hz				X2			
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids			-			К1	-
Light module cold white 6,500 K: light strips arranged on the side walls of the interior, 10 strips for model 110, 14 for model 260/400/750, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)	-				Τ7		-
Light module cold white 6,500 K + warm white 2,700 K: LED light strips - 10 strips for model 110, 14 for model 260/400/750 - (5 resp. 7 alternating cold white light strips and 5 resp. 7 warm white light strips) on the side walls of the interior, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)	-				T8		-
Light module warm white 2,700 K: light strips arranged on the side walls of the interior, 10 strips for model 110, 14 for model 260/400/750, programme-controlled dimming from 0 to 100 % (in 1 % steps), ramp programming in combination with temperature (only with TwinDISPLAY)	-				Т9		-
Interior socket, ampacity 230 V/2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68				R3			
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap, standard positions (F0 and F2 not for model size 260 with light module; F0 - F3 not for model size 110 with light module)left centre/centre right centre/top right centre/top				F0 F1 F2 F3			
Entry port, 23 mm clear diameter, can left be closed by flap (please, state right location) rear				F4 F5 F6			
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)				D6			
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)				F7			
4 - 20 mA current loop interface (-10 to +80 °C = 4 - 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY) - price per sensor				V3 V6			
Works calibration certificate for 3 temperatures: +5 °C, +37 °C, +60 °C				D00129			
Works calibration certificate for one (freely selectable) temperature value according to customer specification				D00109			
Door with lock and key (safety lock)				B6			
Door hinged on the left Potential-free contact (24 V/2 A) with socket, according to NAMUR NE 28 for			B8			-	
external monitoring (indicates when setpoint is reached) Potential-free contact for combination error message (e.g. supply failure, sensor				H5			
fault, fuse)				H6			
Potential-free contact (24 V/2 A) with 2 contacts socket to NAMUR NE 28, for signal generation, controlled by programme segment, for free-selectable functions to be activated (e.g. activation of audible and visual signals, exhaust motors, fans, stirrers, etc.). Only for units with TwinDISPLAY				H72			
Process-dependent programmable door lock (only for units with TwinDISPLAY)				D4			
Door-open-recognition (only for units with TwinDISPLAY) Flexible Pt100 for positioning in chamber or in load with socket, 4-pin, according to NAMUR NE 28, for external temperature recording (load temperature) max. 3 sensors				V5 H4			
Flexible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software				H8			
MobileALERT, notification by SMS in case of any error or alarm of the device. Requires option ${\rm H6}$				C3			
Castor frame (2-part), height 140 mm		R9				-	
Accessories		30	55	110	260 40	00 750	1060
Stainless steel grid, electropolished		E28884	E20164	E20165	E28891	E20182	B41251
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; size 7 and fixing screws (only in connection with option K1). Please consider max. loading	50 with guide bars		-	E29767	E29766	B32190	B32550
Perforated stainless steel shelf		B29727	B03916	B00325	B29725	B00328	B32549
Additional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and in connection with option K1). Please consider max. loading of chamber	fixing screws (only			-		B32191	-

Accessories	30	55	110	260	400	750	1060
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	E02070	E02070 E02072 E02073 E29726			E02075	B32599	
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)			-			B32763	-
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1	B04356	B04358	B04359	B2S	9722	B04362	B29769
Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)			-			B34055	-
Guarantee extension by 1 year		GA1Q5			GA2Q5		GA4Q5
USB-Ethernet adapter				E06192			
Ethernet connection cable 5 m for computer interface				E06189			
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)				B33170			
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number	B33172						
Set of height adjustable feet (4 pcs)		B29	9768			-	
Stacking set (4 pcs) for stacking of appliances of same size		B29744				-	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29734	B29738	B42116	B29	742
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729	B29731	B29735	B29739	B42117	B29	743
Subframe, adjustable in height (size 30 to 75: height 600 mm, size 110 to 450: height 500 mm)	B29745	B29747	B29749	B29751		-	
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746	B29748	B29750			-	
Subframe, adjustable in height, height 130 mm, for example for units with fresh air filter	B33657	B33659	B33661	B33664		-	
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1						
Integration of one additional unit (up to max. 15 units) into an already existent FDA-software licence (only for units with TwinDISPLAY)				FDAQ2			
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer				D00124			
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 - 1060) to DIN 12880:2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer site on demand (GER, AT, CH only)	D00125 D00127						
External measuring instrument with sensors for daylight and UV-light. Product information on demand (models IPPplus)			B04	713			-
External measuring instrument with additional measuring head for temperature and humidity measurement. Product information on demand (models IPPplus)			B04	714			-



Cooled storage incubator IPS with SingleDISPLAY AtmoCONTROL software

Model sizes: 260 / 750 +14 °C to +45 °C

**COOLED STORAGE INCUBATOR IPS** Save energy and reduce the strain on the climate at the same time! If microbiolo-gical cultures, BOB5 samples, drinks containers or cosmetics need to be stored over a long period at constant temperatures, cooled storage incubators IPS with energy-efficient Peltier technology are the perfect choice: absolute reliability, precision, durability and eco-friendliness.





# Considerable potential for savings in acquisition and operating costs

Temperature changes are not always necessary for long-term storage or incubating. So why design heating, cooling and controlling systems for rapid heating up and cooling down times? The performance of the IPS was tailor-made for permanent operation at constant temperatures close to room temperature. The advantage: Acquisition costs and operating costs are considerably reduced in comparison to conventional cooled incubators with compressor technology, as well as to a large Peltier-cooled incubator.

# Ideal for high ambient temperatures

Thanks to Peltier elements integrated for cooling the working chamber, the chamber load won't break into sweat even at high ambient temperatures. Constant and precise incubation at room temperature is guaranteed.

# Low in vibration and durable for absolutely safe long-term storage

Like the cooled incubator IPP, the IPS offers all the advantages of Peltier technology to the user. Its interior chamber that is completely insulated from the environment minimises the risk of drying out of the samples. It is practically noise-free and not only reduces stress on the chamber load but also soothes the nerves of employees thanks to its quiet operation.





Glimpse into a Memmert storage incubator: Peltier elements guarantee perfect climate inside the chamber.

## **COOLED STORAGE INCUBATORS IPS**

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

Interior:	Stainless steel, material 1.4301 (ASTM 304), deep-drawn
Housing:	Textured stainless steel, rear zinc-plated steel, intuitively operated SingleDISPLAY (TFT colour display) with touchscreen
Double doors:	Outside stainless steel, fully insulated, inside glass (size 750 two leaves)
Connection:	Mains cable with plug (German type)
Installation:	4 feet; size 750 mounted on lockable castors
Interfaces:	Ethernet





Standard units are safety-approved and bear the test marks:  $\mathsf{C} \in \mathsf{E}\mathsf{F}\mathsf{I}$ 

Model sizes/Deso	ription		260	750
Stainless steel	Volume	approx. I	256	749
interior	Width	(A) mm	640	1040
	Height	(B) mm	800	1200
	Depth (less 10 mm for fan – Peltier)	(C) mm	500	600
	Max. number of grids/shelves	number	9	14
	Max. loading per grid/shelf	kg	20	30
	Max. loading of chamber	kg	2	00
	Max. loading per slide-in drip tray	kg	4	8
	Max. loading per bottom drip tray	kg	4	8
Textured stainless	Width	(D) mm	824	1224
steel exterior	Height (size 750 with castors)	(E) mm	1183	1726
	Depth (without door handle), door handle + 56 mm	(F) mm	754	856
Standard	Stainless steel grids, electropolished	number		2
equipment	Standard works calibration certificate (measuring point chamber center)	°C	+18 a	nd +25
Temperature	Working temperature range	°C	+14 t	0 +45
	Setting temperature range	°C	+14 t	0 +45
	Setting accuracy	°C	0	.1
Further data	Electrical load at 230/115 V, 50/60 Hz	approx. W	5	50
	Peltier elements in the rear	number		2
Packing data	Net weight	approx. kg	113	230
-	Gross weight (packed in carton)	approx. kg	164	301
	Width	approx. mm	930	1330
	Height	approx. mm	1380	1910
	Depth	approx. mm	930	1050
Order No. Cooled	Storage Incubators			

Order No. Cooled Storage Incubators

IPS260 IPS750

Options	260	7	50
Voltage 115 V, 50/60 Hz		X2	
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids	-	К	.1
Interior socket, ampacity 230 V/2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68		R3	
Entry port, 23 mm clear diameter, for introducing connections at left centre/centre		FO	
the side, can be closed by flap, standard positions left centre/top		F1	
right centre/centre		F2	
right centre/top		F3	
intry port, 23 mm clear diameter, can be closed by flap, in special left		F4	
positions (please state location) right		F5	
rear		F6	
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)		D6	
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please, state location)		F7	
4 - 20 mA current loop interface (0 to +70 °C = 4 - 20 mA)       Temperature controller, actual value		V3	
Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring		V6	
Vorks calibration certificate for one (freely selectable) temperature value according to customer specification	DC	00109	
Door with lock and key (safety lock)		B6	
Door hinged on the left	B8		-
Potential-free contact (24 V/2 A) with socket, according to NAMUR NE 28 for external monitoring (indicates when setpoint is reached)		H5	
otential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)		H6	
lexible Pt100 for positioning in chamber or in load with socket, 4-pin, according to NAMUR NE 28, for external temperature recording load temperature) max. 3 sensors		H4	
Exible Pt100 temperature sensor, positioned flexibly in chamber or load, for local temperature measurement (up to 3 additional sensors are possible). The measured temperature can, if required, be indicated on the display, recorded in the integral data store, and can be documented via the AtmoCONTROL software		H8	
NobileALERT, notification by SMS in case of any error or alarm of the device. Requires option H6		C3	
astor frame (2-part), height 140 mm	R9		-
Accessories		260	750
tainless steel grid, electropolished		E28891	E2018
Additional reinforced stainless steel grid, electropolished, max. loading 60 kg; size 750 with guide bars and fixing screws (only in connection lease consider max. loading of chamber	n with option K1).	E29766	B3219
erforated stainless steel shelf		B29725	B0032
dditional reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (only in connection with option K1). Please bading of chamber	consider max.	-	B3219
tainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1		E29726	E0207
tainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (can be used only in connection with option K1)		-	B3276
tainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution) - cannot be used in connection with option K1		B29722	B0436

Stainless steel bottom drip tray, 15 mm rim (can be used only in connection with option K1)	-	B34055
Guarantee extension by 1 year	GA	2Q5
USB-Ethernet adapter	EOG	5192
Ethernet connection cable 5 m for computer interface	EOG	5189
USB stick with documentation software AtmoCONTROL and operation manual. When reordering please specify serial number	B33	3172
Set of height adjustable feet (4 pcs)	B29768	-
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29738	B29742
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29739	B29743
Subframe, adjustable in height (height 500 mm)	B29751	-
Subframe, adjustable in height, height 130 mm, for example for units with fresh air filter	B33664	-
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer	D00	0124
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 27 measuring points to DIN 12880;2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer	D0(	0127

measuring points to DIN 12880:2007-05. PQ check list as support for validation by customer. Price for further temperature values and validation at customer site on demand (GER, AT, CH only)

## **MODEL VARIANTS**

SingleDISPLAY ControlCOCKPIT with one TFT display	TwinDISPLAY ControlCOCKPIT with two TFT displays					
AVAILABLE APPLIANCES UN/UNm / UF/UFm / IN/INm / IF/IFm / IFbw / SN / SF / IPP / IPS	AVAILABLE APPLIANCES UNplus/UNmplus / UFplus/UFmplus / UF TS / UNpa INplus/INmplus / IFplus/IFmplus / SNplus / SFplus / VO ICOmed / IPPplus / ICPeco / ICP / HPP / ICHeco / ICH / HCP					
One high-resolution TFT colour display with touch-sensitive buttons for selection of functions	Two high-resolution TFT colour displays with touch-sensitive buttons for selection of functions					
Available parameters on the ControlCOCKPIT: Temperature (Celsius or Fahrenheit), fan speed, exhaust air flap position, programme time	Available parameters on the ControlCOCKPIT: Temperature (Celsius or Fahrenheit), fan speed, exhaust air flap position, programme time, relative humidity, illumination, CO <sub>2</sub>					
One temperature sensor Pt100 DIN class A in a 4-wire circuit	Two Pt100 sensors DIN class A in a 4-wire circuit for mutual monitoring, taking over functions in case of an error					
	HeatBALANCE function for application specific adjustment of heat output distribution (balance) between the upper and lower heating groups in an adjustment range between -50 % and +50 % (not valid for models 30, HPP110, IPP110plus, ICP, ICH)					
AtmoCONTROL software for reading out, managing and organising the data logger via Ethernet interface (temporary trial version can be downloaded). USB stick with AtmoCONTROL software available as accessory (on demand)	AtmoCONTROL software on a USB stick for programming, managing and transferring programmes via Ethernet interface or USB port					
	ControlCOCKPIT with USB port for uploading programmes, reading out protocol logs, activating the User-ID function					
	Displaying of already logged protocol data on the ControlCOCKPIT (max 10,000 values correspond to approx. 1 week)					
Ethernet interface on the rear of the appliance for reading out the protocol log and for online logging	Ethernet interface on the rear of the appliance for reading out the protocol log and for uploading programmes and for online logging					
Double overtemperature protection: Electronic temperature monitoring with freely adjustable monitoring temperature, for models U, I, S with option A6 TWW/TWB (protection class 3.1 or 2), mechanical temperature limiter TB acc. to DIN 12880	Multiple overtemperature protection: Electronic temperature monitoring TWW/TWB (protection class 3.1 or 2 resp. 3.3 for units with active cooling) and mechanical temperature limiter TB (protection class 1) acc. to DIN 12880, AutoSAFETY automatically adjusts to the set value within a freely adjustable tolerance range. Setting individual MIN / MAX values for over/undertemperature and also for all other parameters such as relative humidity, CO <sub>2</sub>					
PID microprocessor control with	integrated auto-diagnostic system					
Structured stainless steel housing, scratch-resist	tant, robust and durable; rear of zinc-plated steel					
	ar of the appliance for single-phase power specific systems and IEC standards					
Internal data logger with a stor	rage capacity of at least 10 years					
German, English, French, Spanish, Polish, Czech, Hungarian language settings available on the ControlCOCKPIT						
Digital backwards counter with target time	setting, adjustable from 1 minute to 99 days					

The SetpointWAIT function guarantees that the process time does not start until the set temperature is reached at all measuring points – optional for temperature values recorded by the freely positionable Pt100 sensors inside the chamber

Adjustment of three calibration values for temperature and additional appliance specific parameters directly at the ControlCOCKPIT

## SOFTWARE AtmoCONTROL


## AtmoCONTROL

## The innovative control and logging software

Parameters such as temperature and humidity as well as the process time can be set directly at the ControlCOCKPIT. Ramp programming is done via the control and logging software AtmoCONTROL, which features a completely new software design.

## Drag, drop & go!

Numerical and graphic programming of complex processes is a thing of the past. Today, programming is done via AtmoCONTROL by means of the mouse or touchpad on your notebook. Even the most complex ramp programmes are created within minutes. Simply drag & drop the graphical symbols for the desired parameters to the input field and change the values according to your wishes with a mouse click.




# Programme functions for appliances with SingleDISPLAY and TwinDISPLAY

- · Reading out, managing and organising the data logger
- Saving the log memory in various formats
- Online monitoring of up to 32 connected appliances
- Optical alarms when the alarm limits individually set at the ControlCOCKPIT are exceeded
- Automatic alarm to one or several e-mail addresses

## Additional functions for appliances with TwinDISPLAY

- Intuitive programming and archiving of ramps and programme sequences
- Synchronous visualisation of the created programme sequence during programming
- Application-specific repeat functions (loops) can be inserted within a temperature control programme in any place
- Simple creation of repeating weekly programmes
- Programming, managing and transferring programmes via Ethernet interface or USB port





## HEATING AND DRYING OVENS

- UNIVERSAL OVEN U
- PASS-THROUGH OVEN UF TS
- PARAFFIN OVEN UNpa
- STERILISER S
- VACUUM OVEN VO
- BLANKET WARMER IFbw

## **INCUBATORS**

- INCUBATOR I
- CO, INCUBATOR ICOmed
- COMPRESSOR-COOLED INCUBATOR ICPeco/ICP
- PELTIER-COOLED INCUBATOR IPP
- COOLED STORAGE INCUBATOR IPS

## CLIMATE CHAMBERS

- CONSTANT CLIMATE CHAMBER HPP
- HUMIDITY CHAMBER HC
- CLIMATE CHAMBER ICHeco/ICH
- ENVIRONMENTAL TEST CHAMBER CTC/TTC

## WATERBATHS / OILBATHS

- WATERBATH W
- OILBATH O



Memmert GmbH + Co. KG P.O. Box 1720 | D-91107 Schwabach Tel. +49 9122 925-0 | Fax +49 9122 14585 E-Mail: sales@memmert.com facebook.com/memmert.family The platform for experts: www.atmosafe.net