

# GVPC AGAR FOR *LEGIONELLA*

DETECTION AND ENUMERATION OF *LEGIONELLA*

BM07108

## 1 INTENDED USE

GVPC agar for *Legionella* is used for the enumeration, isolation and culture of *Legionella* species in water and other samples susceptible of harboring the bacteria.

## 2 HISTORY

In 1977, MacDade *et al.* Were the first to isolate the agent responsible for Legionnaire's Disease, a bacteria now known as *Legionella pneumophila*. After this discovery, numerous occurrences of *Legionella* isolation of *Legionella* were reported in fresh water environments such as water distribution systems, air conditioning, cooling towers, and spas. 48 species of *Legionella* are currently known.

In 1978, Weaver succeeded in cultivating *Legionella* on Mueller-Hinton chocolate agar. Feeley *et al.*, deduced that cysteine and ferric pyrophosphate could replace the vitamin and hemoglobin supplements found in the Mueller Hinton chocolate agar. Their work led to the formulation of a medium dubbed F-G agar. They determined as well that an atmosphere enriched at 2.5% CO<sub>2</sub> was necessary for *Legionella* culture. In 1979, Feeley *et al.* modified the the F-G medium by replacing acid hydrolysate of casein by yeast extract, and adding activated charcoal while eliminating starch. The resulting CYE media allowed better growth of *Legionella*. In 1980, Pasculle *et al.* supplemented the CYE medium with ACES buffer. They demonstrated that this new medium, designated BCYE, offered a better recovery of *Legionella* and could be incubated aerobically. In 1981, Edelstein increased the sensitivity of the medium by adding  $\alpha$ -cetoglutarate (BCYE $\alpha$  medium), and Wadowsky & Yee suggested incorporating glycine, vancomycin and polymyxin B (GVP medium) to obtain a selective culture media. In 1984, Dennis *et al.* formulated the current GVPC medium by adding cycloheximide into GVP medium. They demonstrated that this selective medium allowed a greater level of *Legionella* isolation.

## 3 PRINCIPLES

Yeast extract constitutes a primary nutrient leading to *Legionella* growth.

Activated charcoal decomposes hydrogen peroxide (toxic metabolic by-product), captures the carbon dioxide and modifies the surface tension.

The ACES/KOH buffer maintains the pH and permits aerobic incubation.

Cysteine and ferric pyrophosphate represent indispensable nutritive elements for the growth of *Legionella*.

$\alpha$ -cetoglutarate is a growth activator for *Legionella*.

Secondary microflora are inhibited by the association of glycine, vancomycin , polymyxin B and cycloheximide.

## 4 INSTRUCTIONS FOR USE

### Surface inoculation :

- To the surface of pre-poured plates (BM071), or to plates prepared as above, transfer 0.1 mL of the sample to be tested and its serial dilutions.
- Spread the inoculum with a sterile triangle or "hockey stick".
- Incubate at 36 ± 2 °C for 10 days.
- During the incubation period, examine the plates starting from the third day and three successive moments at intervals of 2 to 4 days.

✓ **Inoculation :**  
**On surface**  
**After membrane filtration**

✓ **Incubation :**  
**10 days at 36 ± 2°C**

### After membrane filtration :

- Aseptically filter through a membrane a known volume of the sample to test (10 mL to 1000 mL of water sample).

- Deposit the membrane on the surface of pre-poured plates (BM071), filtered side up and making sure that the membrane and agar are in close contact. The plates should be brought to room temperature before use.
- Incubate at  $36 \pm 2$  °C for 10 days.
- During the incubation period, examine the plates starting from the third day and three successive moments at intervals of 2 to 4 days.

## 5 RESULTS

Colonies of *Legionella* spp. present a white to gray coloration. They can also have blue, pink, purple, maroon, greenish-yellow or dark red pigmentation that fades, becoming whiter and filamentous with age. Their surface is smooth with precise edges. Some strains may give a ground glass or “fried egg” aspect when observed through a binocular scope, while others may present a brilliant white fluorescence under a UV light.

Colonies of *Legionella* that develop on white membrane filters may have a different aspect to those that develop against a black background.

Enumerate each colony type separately. Select at least three characteristic colonies of *Legionella* on each of the agar plates. Re-streak each colony onto a plate of BCYE $\alpha$  without cysteine (BM073) and a plate of BCYE $\alpha$  (BM072) from the BT007 kit.

See ANNEX 1 : PHOTO SUPPORT.

## 6 TYPICAL COMPOSITION

The typical composition can be adjusted to obtain optimal performance.

For 1 liter of medium :

- Yeast extract .....	10.0 g
- Activated charcoal.....	2.0 g
- $\alpha$ -cetoglutamate, monopotassium salt.....	1.0 g
- ACES (2-[2-amino-2-oxoethyl)-amino] ethanesulfonic acid).....	10.0 g
- Potassium hydroxide.....	2.8 g
- L-cysteine, hydrochloride .....	0.4 g
- Ferric pyrophosphate .....	0.25 g
- Glycine .....	3.0 g
- Vancomycin.....	1.0 mg
- Polymyxin B.....	80000 IU
- Cycloheximide.....	80.0 mg
- Bacteriological agar.....	12.0 g

pH of ready-to-use medium at 25°C :  $6.9 \pm 0.1$ .

## 7 QUALITY CONTROL

Prepared media in plates : black agar, with visible particles of activated charcoal.

Typical cultural response after 5 days incubation at 36 °C (NF T 90-461) :

Microorganisms		Growth
<i>Legionella pneumophila</i>	WDCM 00107	$66\% \leq R_2 \leq 150\%$
<i>Enterococcus faecalis</i>	WDCM 00176	Inhibited
<i>Escherichia coli</i>	WDCM 00179	Inhibited

## 8 STORAGE

**Pre-poured media in Petri dishes** : 2-8°C, shielded from light.

Expiration dates are indicated on the labels.

## 9 PRESENTATION

**Pre-poured Petri dishes (Ø 90 mm)**

20 plates ..... BM07108

## 10 BIBLIOGRAPHY

---

McDADE, J.E., SHEPARD, C.C., FASER, D.W. et al.. 1977. Legionnaires' disease : isolation of a bacterium and demonstration of its role in other respiratory disease. *New England Journal of Medicine*, **297** : 1197-1203.

FEELEY, J.C., GORMAN, G.W., WEAVER, R.E., MACKEL, D.C. and SMITH, H.W.. 1978. Primary isolation media for Legionnaires disease bacterium. *Journal of Clinical Microbiology*, **8** : 320-325.

WEAVER, R.E. and FEELEY, J.C. 1979. Cultural and biochemical characterization of Legionnaires' disease bacterium. In *Legionnaires', the Disease, the Bacterium, and Methodology*. Centers for Disease Control, Atlanta, GA. pp 20-25.

FEELEY, J.C., GIBSON, R.J., GORMAN, G.W., LANGFORD, N.C., RASHEED, J.K., MACKEL, D.C. and BAINE, W.B.. 1979. Charcoal-yeast extract agar: primary isolation medium for *Legionella pneumophila*. *Journal of Clinical Microbiology*, **10** : 437-441.

PASCULLE, A.W., FEELEY, J.C., GIBSON, R.J., CORDES, L.G., MYEROWITZ, R.L., PATTON, C.M., GORMAN, G.W., CARMACK, C.L., EZZELL, J.W. and DOWLING, J.N.. 1980. Pittsburgh pneumonia agent: direct isolation from human lung tissue. *Journal of Infectious Diseases*, **141** : 727-732.

EDELSTEIN, P.H.. 1981. Improved semiselective medium for isolation of *Legionella pneumophila* from contaminated clinical and environmental specimens. *Journal of Clinical Microbiology*, **14** : 298-303.

WADOWSKY, R.M. and YEE, R.B.. 1981. Glycine-containing selective medium for isolation of *Legionellaceae* from environmental specimens. *Applied and Environmental Microbiology*, **42** : 768-772.

DENNIS, P.J.L., BARTLETT, C.L.R. and WRIGHT, A.E.. 1984. Comparison of isolation methods for *Legionella* spp. In Thornsby, C. *et al.* (ed.) *Legionella : Proceedings of the 2<sup>nd</sup> International symposium Washington D.C.* American Society for Microbiology, p 294-296.

Journal Officiel du 31 décembre 2004. Arrêté du 13 décembre 2004 relatif aux installations de refroidissement par dispersion d'eau dans un flux d'air soumises à autorisation au titre de la rubrique n° 2921.

NF T 90-431 September 2003 and NF T 90-431/A1 April 2006. Water quality. Detection and enumeration of *Legionella* spp and *Legionella pneumophila*. Method by direct inoculation and after concentration by membrane filtration of centrifugation.

FD T 90-522. July 2006. Water quality. Sampling technical guide for the detection of *Legionella* in water.

NF EN ISO 11731-2 (T 90-430). July 2008. Water quality. Detection and enumeration of *Legionella*. Part 2 : direct membrane filtration method for waters with low bacterial counts.

## 11 ADDITIONAL INFORMATIONS

---

The information provided on the package take precedence over the formulations or instructions described in this document and are susceptible to modification at any time, without warning.

Document code : BM071/EN/2004-06 : 5.  
Creation date : 2004-06  
Update : 2013-06  
Grounds for revision : General revision (§ 7 : Quality control).

## GVPC Agar for *Legionella*

Detection and enumeration of *Legionella*.

### Reading :

Growth obtained after 10 days of incubation at 36 °C.



### *Legionella pneumophila*

Characteristic colony  
White to gray color with a smooth surface ;  
some colonies may present a ground glass  
appearance under a binocular scope