

It does the counting, so you can do the work that counts.

Introducing the Neogen® Petrifilm® Automated Feeder.

Enumeration is a critical but time-consuming task. And you've got your best people on it because they're the most reliable. But what if you could free up almost half their day? It could take a technician up to 4 hours to count 300 plates. The Petrifilm Automated Feeder does it in 33 minutes. Ideal for labs that process more than 100,000 Petrifilm Plates a year, the Petrifilm Automated Feeder increases capacity and productivity, allowing you to meet more customer needs with less labor.

Meet your new assistant.

It's always reliable, it decreases technician-to-technician variability, and it doesn't need training. And we must say, it comes with an impressive résumé.

- Automatically feeds and enumerates up to 300 Petrifilm Plates in any order in 33 minutes.
- Built to be robust, the Petrifilm Automated Feeder has a separation station that picks and separates plates, ensuring continuous throughput, and your peace of mind.
- Continuously pushes data into an Excel file for seamless LIMS integration. No more transcription!



Learn how the Petrifilm Automated Feeder can increase productivity.

info.neogen.com/PetrifilmAutomation

More reliability. Less variability.

Compatibility

Same easy-to-use software as Petrifilm® Plate Reader Advanced with added features for Petrifilm Automated Feeder.


Works with all Petrifilm Plate Reader Advanced devices.

Compatible with 10 types of barcodes; plates can be read in any order.



Separation Station

Separates plates that stick together.




Sorting Station

Groups plates based on customizable preferences.



Automated Enumeration

The input bin holds 300 Petrifilm Plates, which can be enumerated in 33 minutes.



Automated Transcription

Automatically pushes data into a Excel file for seamless LIMS integration.

Results can now be reviewed in a thumbnail view, which allows for easy screening of hundreds of tests at a time.

Ordering Information

SKU	Description
700006589	Petrifilm Automated Feeder, 110V
700006590	Petrifilm Automated Feeder, 220V
700006593	Petrifilm Automated Feeder, Input Bin
700006594	Petrifilm Automated Feeder, Output Bin

SKU	Description
700006587	Petrifilm Automated Feeder, Station Cups, 5-pk
700006588	Petrifilm Automated Feeder, Picker Cups, 5-pk
700006592	Petrifilm Automated Feeder, Filters, 4-pk
700006591	Petrifilm Automated Feeder, Extended Service Agreement

