

BagSeal®

Thermal sealing unit

Ref 261 000



Manufactured under:



BagSeal® allows a wide and tight sealing of all kinds of bags (with/without filter). Perfect closure is guaranteed.

TECHNICAL SPECIFICATIONS

- All stainless steel except handles, welding blades (anodized aluminum), bar joint (chromated galvanized steel)
- Adjustable power heating
- LED light code at the end of sealing (welding process duration indicator: welding + cooling time)
- Large sealing: 250 x 5 mm
- For bags 80/100/400/2000 mL
- Pressure: 22 Newton/cm²

- Voltage and frequency: 100-127 V ~ / 220-240 V ~ 50-60 Hz
- Max power: 1600 W

- International protection marking: IP40
- Dimensions (w x d x h): 35.9 x 27.6 x 32.2 cm, weight: 10.8 kg
- Box (w x d x h): 49 x 40 x 49 cm, gross weight: 13.6 kg

- Warranty: 1 year further to warranty card recording
- Availability of spare parts: 10 years
- Manufactured under CE, UKCA, RoHS, WEEE
- Designed and made in France

OPTIONS / ACCESSORIES

- Ref 261 500 Maintenance kit for BagSeal®: 2 sealing threads + 2 fuses + 2 Teflon strips