

## BagClip®

### Closing clip

Ref 231 010 (80/100 mL) - Ref 231 040 (400 mL) - Ref 231 350 (3500 mL)



Manufactured under:



BagClip® allows a temporary closing in 3 points of any type of bag, except for the Rollbag range.

### TECHNICAL SPECIFICATIONS

- Gas and liquid leakproof
- No contact with the sample
- Available for 80/100, 400, 2000/3500 mL blender bags
- Designed and made in France

#### BagClip® 80/100

Ref: 231 010

- Available for 80- and 100-mL bags
- Length (white clip): 130 mm
- Length (green tube): 160 mm
- Net weight: 0.009 kg
- Box (w x d x h): 21.5 x 16 x 6.5 cm, gross weight: 0.74 kg
- Quantity per pack: 50

#### BagClip® 400

Ref: 231 040

- Available for 400 mL bags
- Length (white clip): 195 mm
- Length (green tube): 225 mm
- Net weight: 0.013 kg
- Box (w x d x h): 28.5 x 19.5 x 6.5 cm, gross weight: 0.77 kg
- Quantity per pack: 50

#### BagClip® 3500

Ref: 231 350

- Available for 2000 and 3500 mL bags
- Length (white clip): 400 mm
- Length (green tube): 450 mm
- Net weight: 0.027 kg
- Box (w x d x h): 50.5 x 11.5 x 10.5 cm, gross weight: 1.57 kg
- Quantity per pack: 50

### CAUTION

- Do not put in an autoclave.
- Do not put in a washing machine.
- Can be disinfected with ethanol 70 %, hydrogen peroxyde 6 % or bleach 10 % max.