

DAILY SHIFT HEAD

TRIO.BAS™

ORUM INTERNATIONAL
© ALL RIGHTS RESERVED

Sterile aspirating head for TRIO.BAS air samplers



340 (certification of sterilization in each box)



341 (certification of sterilization in each box)



DESCRIPTION

- The sterile Daily Shift aspirating heads (DSH - sterile Daily Shift Head) avoid the sterilization process necessary for stainless steel aspirating heads.
- The sterilization is proven by an official certificate. This document is requested by regulatory authorities.
- The double irradiated sterile packaging allows the users to always have aspirating heads ready for use.
- The transparency of the sterile Daily Shift is useful for verifying the culture plate is inserted correctly in the aspirating chamber.
- Main customers for sterile Daily Shift heads are: agro-food industries, dairy, catering, HACCP, beverage, cosmetic, sewage treatment plants, outdoor environments, primary and secondary schools, pharmaceuticals, cleanrooms, biotech, hospitals, clinics, microbiological labs, HVAC building monitoring, environmental labs, healthcare ambient monitoring, health authorities.
- They are suitable for all TRIO.BAS air samplers.
- Shelf life: 6 years from the date of sterilization.

IDENTIFICATION CODES

Code	Sterile Aspirating head
340	Sterile Daily Shift Aspirating Head for Contact plate 55 mm - in double sterile bag (27xbox)
341	Sterile Daily Shift Aspirating Head for Petri 90 mm plate - in double sterile bag (27xbox)

