



FLUICS CRYO LABEL SPECS

Smear proof	Spraying and swabbing with alcohol and surface sanitizing wipes (ethanol, isopropanol, DMSO, bleach)
Adhesive	Extra permanent, especially on curved surfaces (Eppis, cryo tubes); adhere to wet and moist surfaces
Temperature range	-269°C to +110°C incl. liquid He and N2
Sterilization	Ethylene gas sterilization and gamma-irradiation
High contrast readability	Black text on yellow labels
Size	25 mm x 13 mm (1 inch x 0.5 inch)
Easy handling with gloves	✓
Compliant for indirect food contact	FDA 21 CFR 175.105
Printing supplies (ribbons)	Included in label price