


## 50mL Polypropylene Conical Centrifuge Tubes with Printed Graduations

Item# 6288/6289/6287



Item# 6288

### Description

General purpose centrifuge tube, 500/case. 50mL capacity. Blue screwcap. Black printed graduations. White printed writing space. Available in sterile with attached caps and non-sterile with separate caps.

### Features

- Tube molded from chemical resistant, hydrophobic polypropylene (PP)
- Cap molded from polypropylene (PP)
- Conical style bottom
- Printed black graduations (0.5mL graduations)
- White writing area

### Products

Item#	Description	Packaging
<b>6288</b>	50mL, PP, sterile, 25 per bag	500 per case
<b>6289</b>	50mL, PP, sterile, 25 per rack	500 per case
<b>6287</b>	50mL, PP, non-sterile, bulk, with separate blue screwcap	500 per case



Item# 6289 in 25 place rack

## 50mL Polypropylene Conical Centrifuge Tubes with Printed Graduations


Item# 6288/6289/6287

### Technical Specifications

Material: Tube: Hydrophobic Polypropylene (PP)  
Cap: Polypropylene (PP)  
Sterilization: Ethylene Oxide (EO)  
Autoclavable: Yes (without cap)  
Temperature range: -20°C to +121°C  
Spin rate: 8500 RCF

### Product Certifications

- ISO 13485
- CE



#### Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail [mail@globescientific.com](mailto:mail@globescientific.com)

#### Globe Scientific Inc.

610 Winters Avenue • Paramus, NJ 07652

[www.globescientific.com](http://www.globescientific.com)

© Copyright 2021 Globe Scientific Inc. All rights reserved.