

PRODUCT INFORMATION

Beef Extract

Cat. No. B02-104

DESCRIPTION

Beef Extract is used as a nutrient source in various culture media. It is prepared from very low-fat bovine skeletal muscle and marrow.

PHYSICAL AND CHEMICAL CHARACTERISTICS

TEST TYPE	SPECIFICATION
Appearance	Beige powder
Total Nitrogen	≥ 10.0%
Loss on Drying	≤ 6.0%
Ash	≤ 15.0%
pH (2% solution)	6.5 – 7.5

STORAGE

Store the sealed bottle containing the dehydrated medium at 2 to 30°C. Once opened and recapped, place the container in a low humidity environment at the same storage temperature. Protect it from moisture and light. The dehydrated medium should be discarded if it is not free flowing or if the color has changed from the original color.