

BagMixer® 400 CC

400 mL lab blender

Ref 024 230

WINDOW DOOR
LIFETIME
WARRANTY

SHOCK ABSORBERS
LIFETIME
WARRANTY

3
YEAR
WARRANTY



Standards:



Manufactured under:



Product made for INTERSCIENCE by Interlab, an ISO 9001 certified company

BagMixer® 400 CC is an easy-to-use and powerful lab blender. Adapted to all kinds of applications and with a guarantee of optimal bacterial extraction, BagMixer® 400 CC, equipped with its window door and patented adjustable paddles, ensures fast and perfect homogenization for lab analyses. Cleaning is easy with Click and Clean® removable paddles system.

TECHNICAL SPECIFICATIONS

- All stainless steel 304L
- Window door
- Multi-function digital screen
- Optimal blending volume: 50 - 400 mL
- Variable speed (4 - 6 - 8 - 10 strokes/s)
- Variable blending time (1 s – 1 h or ∞ / count down)
- Adjustable blending power (patented)
- Click & Clean® removable paddles (patented)
- Side-by-side paddle stop
- Q-Tight® closing: tight closing for no more leaks from the top of the bag
- Fully opening door (270°)
- Silent brushless motor
- Security drip tray

- Voltage and frequency: 100-240V~ 50-60Hz
- Fuses: T3.15A L 250V 5x20 mm
- Max power: 55 W

- Dimensions (w x d x h): 26.5 x 42 x 27.5 cm, weight: 16 kg
- Box (w x d x h): 49.5 x 35 x 36.5 cm, weight: 18.6 kg

- 3-year warranty
- Shock absorbers lifetime warranty
- Window door lifetime warranty

- In compliance with: ISO 7218, ISO 6887-1 and FDA BAM (Bacteriological Analytical Manual)
- Manufactured under CE, RoHS, WEEE - Products made for INTERSCIENCE by Interlab, an ISO 9001 certified company
- Designed and made in France

DELIVERED WITH

- 1 security drip tray
- 1 power cable
- 1 user's manual
- 1 quick user guide
- 1 free pack of blender bags